

TECHNICAL DATA

TRIAL FRAMES

	ESS 304	ESS 204
CAPACITY	10 lenses	10 lenses
PUPILLARY DISTANCE	46-82 mm	48-80 mm
WEIGHT	64 g	66 g
OPTIONAL	Polarizing filters ref: ESS305	-

TRIAL CASES

	ESS 260	ESS 167
SPHERIC LENSES	±0.25 to ±20.00 D	±0.25 to ±16.00 D
CYLINDRIC LENSES	±0.25 to ± 6.00 D	±0.25 to ± 4.00 D
PRISMATIC LENSES	0.5 to 10 D	1 to 5 D
DIMENSIONS	Tray only is L 440 x D 320 x H 45 mm Full box is L 470 x D 350 x H 100 mm	Tray only is L 440 x D 265 x H 45 mm Full box is L 465 x D 300 x H 100 mm
WEIGHT	4,750 kg	3.8 kg

Trial frames Trial cases

Comelink Communication Ref: PID06900E - 09/2014 - In the interest of improvement, these specifications are not legally-binding and may be modified without prior notice.

MAKING REFRACTION SIMPLER AND ACCURATE



ESS RANGE

TRIAL FRAMES AND TRIAL CASES

Your contact:


MAKING REFRACTION SIMPLER AND ACCURATE

A complete range of optometric tools to define with precision the prescription needs of your patients.



ESS 304 ERGONOMIC DESIGN

- Easy to manipulate with flexible, smooth turning dials.
- Adjustable nose bridge and temples (length and angle) to adapt to all facial morphology.
- Special grooves to accommodate up to 10 lenses.
- Temples designed to prevent hair entanglement.
- **Optional:** pair of polarized filters for far and near vision examinations.



Pair of filters
Ref: ESS305

ESS 204 FOR UNIVERSAL USE

- Simple and flexible to adjust when patients wearing trial lenses.
- Adjustable nose bridge and temples (length and angle).
- Grooves that allow testing of up to 10 lenses .



ESS 260 FULL DEMO LENS KIT

- Complete kit with 260 lenses.
- Covers 22mm diameter range.
- Accessories: cross, Maddox, cross cylinders, occluders, red/green filters, plano lenses, stenopeic hole.
- Removable tray for easy placement in all dispensing areas.



ESS 167 DEMO LENS SELECTION

- Kit of 167 lenses with essential accessories for measuring visual acuity.
- Covers 22mm diameter range.
- Accessories: occluders, cross cylinders, slots, red/green filters, stenopeic hole.
- Removable tray.

SELECTION OF OPTOMETRIC ACCESSORIES



ESSILOR STEREO TEST to screen for binocular vision problems.

MAG200



ISHIHARA TESTS to screen for color vision anomalies.

14 sheets	MAG301
24 sheets	MAG302
38 sheets	MAG303



PRISMS FOILS (to 67 mm) to compensate deviations.

Pack of 2 or single unit.